infiTOF Application Note In situ monitoring of reaction on MOVPE by using "infiTOF"



In situ monitoring of Ga(CH₃)₃ reaction with HCl on nitride semiconductor Meal - Organic Vapor Phase Epitaxy (MOVPE) by using infiTOF

This application includes the experimental data provided by Dr. Amano, Dr. Nitta, Dr. Nagamatsu and Mr. Ye of Institute of Materials and Systems for Sustainability, Nagoya University. MSI.TOKYO,Inc.(Kanomax Group) has put together this based on the helpful discussions with them.



Introduction

Nitride semiconductors are candidate materials for high-power transistors. To achieve high breakdown voltage performance, the GaN drift layer must be grown with the lowest amount of impurities possible. To improve device performance, several researchers have focused on reducing the carbon, silicon, and oxygen impurities in GaN. Reducing the carbon concentration is especially difficult under conventional growth conditions, because using trimethyl gallium (TMG), which provides a methyl group needed in the production process, also results in carbon impurities. Metalorganic halide-vapor phase epitaxy (MOHVPE) can effectively reduce carbon incorporation by replacing the carbon-based methyl groups with chlorine using HCl. Previous reports on GaN growth by MOHVPE have focused on high growth rate of the GaN bulk crystal without considering the effects of impurity incorporation. However, Amano et al. reported on the impurity concentration in GaN grown by MOHVPE and concluded that chlorine replacement cannot sufficiently reduce the carbon concentration and that the direct reactions must be monitored in the vapor phase.

Experimental

Figure 1 shows the experimental setup of the MOVPE reactor coupled to the infiTOF. In the figure, all gases flow from left to right and then out the exhaust port. TMG reacts with HCl on the wafer surface. The gas products are sampled by the line installed at the center of the wafer's upper surface, and introduced into the electron ionization (El) source of the infiTOF by a microtube (I.D. 0.1mm, L50mm) heated at 120 °C in order to prevent sample gases from adsorbing to the tube. The size of the quartz flow channel in the reactor was designed for a single 4-inch wafer. The number of turns for the infiTOF was set to achieve a resolving power of 10,000 or more. Other MOVEP experimental conditions are listed in Table 1.

Wafer surface temperature	RT-1150 ℃
Wafer material	Sapphire 2"
Channel width	120 mm
Channel height	6 mm
Reactor pressure	200 torr
Flow speed	1 m/sec
Carrier gas	H 2(+N2)
TMG	100 sccm




Fig. 1 Experimental setup of MOVPE reactor coupled to the infiTOF.



Results and Discussion

Figure 2 shows the mass spectrum of TMG in the absence of HCl supply. Peaks for $Ga(CH_3)_3^+$, $Ga(CH_3)_2^+$, $GaCH_3^+$, Ga^+ , and CH_4^+ were observed. The strongest peak was for $Ga(CH_3)_2^+$ while $Ga(CH_3)_3^+$ could almost not be detected at all. The signal intensity ratio of $Ga(CH_3)_3^+$ / $Ga(CH_3)_2^+$ was less than 1%.

Table 2 shows the theoretical and calculated values of each component and their associated errors, which were obtained by a post-calibration process using HCl and Ga(CH3)2 as calibrants. In each component, the mass error between theoretical and calculated values was within 3 [ppm], indicating that there is sufficient mass accuracy to identify the components by exact mass.

each component and their errors.					
	Theoretical	Experimental	u	ppm	
Ga(CH ₃)Cl ₂	155.8849	155.8847	0.00025	1.6	
GaCl ₂	140.8614	140.8611	0.00032	2.3	
Ga(CH ₃) ₂ Cl	133.9414	133.9414	0.00002	0.1	
Ga(CH ₃)Cl	118.9179	118.9179	0.00004	0.3	
GaCl	103.8944	103.8946	0.00015	1.4	

Table 2 Theoretical and calculated values of each component and their errors.




Fig. 2 Mass spectrum of TMG.




Fig. 3 Mass spectrum of reaction products of TMG and HCl.

Regarding the generation of $Ga(CH_3)xCly$ (x = 0, 1, 2. y = 1, 2.) compounds, there are several possibilities besides generation by the reaction of TMG with HCI:

- GaCl resulted from El fragmentation of Ga(CH3)Cl2, Ga(CH3)2Cl, GaCl2, and Ga(CH)Cl.
- Ga(CH3)Cl resulted from El fragmentation of Ga(CH3)Cl2 and Ga(CH3)2Cl.
- GaCl2 resulted from El fragmentation of Ga(CH3)Cl2.



In order to reveal the origin of Ga(CH3)xCly products, mass chromatograms of each component normalized to HCl signal intensity were measured (Figure 4).

Figure 4 shows that:

(1) The signal intensity of Ga(CH3)2Cl and Ga(CH3)Cl clearly decreased.

- (2) The signal intensity of GaCl clearly increased.
- (3) The signal intensity of Ga(CH3)Cl2 and GaCl2 seem to decrease somewhat.

Because only GaCl increased, it can be concluded tha GaCl was not formed by El fragmentation of Ga(CH₃)Cl₂, Ga(CH₃)₂Cl, Ga(CH₃)₂Cl, or GaCl₂, but by reaction of TMG with HCl. Mass chromatograms of Ga(CH₃)₂Cl and Ga(CH₃)₂Cl had decreasing trends, so it is possible that Ga(CH₃)₂Cl was the result of El fragmentation of Ga(CH₃)₂Cl. It is difficult to discuss the origin of GaCl₂ and Ga(CH₃)₂Cl₂, because the associated mass chromatograms do not exhibit a concrete trend. Since the origin of Ga(CH₃)₂Cly products is based only on experimental results, it is necessary to supplement the results with theoretical calculation in order to precisely understand the reaction of TMG and HCl.




Fig.4 Normalized mass chromatograms of chlorine adducts.

Conclusion

In-situ monitoring of the reaction of TMG with HCl on nitride semiconductor MOVPE was performed using infiTOF. Several reaction products such as GaCl, Ga(CH3)Cl2, Ga(CH3)2Cl, Ga(CH3)Cl, and GaCl2 were observed. Results indicate that GaCl resulted from reaction of TMG with HCl, and Ga(CH3)Cl resulted from El fragmentation of Ga(CH3)2Cl. It is necessary to compare experimental results with theoretical results in order to understand the true nature of the reaction of TMG with HCl in the MOVPE reactor. Results also indicate that infiTOF is useful for investigating MOVPE reactions.

Copyright© 2017 MSI.TOKYO, Inc.



KANOMAX CORPORATION

<Oosaka Head Office>
2-1 Shimizu Suita-shi, Osaka 565-0805 JAPAN
TEL: +81-6-6877-0463 / FAX: +81-6-6879-2080
<Tokyo Office>
2-6-2 Hamamatsu-cho, Minato-ku, Tokyo 105-0013 JAPAN
TEL: +81-3-5733-6544 / FAX: +81-3-5733-6545
E-mail: analytical@kanomax.co.jp
URL: https://kanomax.biz/asia/
AANMT-009-2E