

Product	Model/Part Number	Unit	Shelf Life (months)	Refrig Storage?	Comment
LNS Volume Standard (5E17 nm <sup>3</sup> /mL)	KFMT LNS Vol Std	125 mL bottle	12	N	
Scanning TPC Calibration Kit consisting of:	KFMT STPC Cal Kit	8-30 mL bottles	3	Y	
30 nm silica - 1E11/mL > 10 nm in 30 mL bottle	KFMT CS 30nm 1e11	2-30 mL bottles	6	Y	if purchased separately
30 nm silica - 1E10/mL > 10 nm in 30 mL bottle	KFMT CS 30nm 1e10	2-30 mL bottles	3	Y	"
UHPW in 30 mL bottle	KFMT UHPW	2-60 mL bottles	3	Y	"
Poly (10, 12, 20, 30 nm) silica in 30 mL bottle	KFMT PSL Cocktail 10-30nm	2-30 mL bottles	6	Y	"
20-125 nm Poly PSL (5E8/mL > 30 nm nominal)	KFMT cocktail 20-125nm	125 mL bottle	3	Y	
20-220 nm Poly PSL (5E8/mL > 30 nm nominal)	KFMT cocktail 20-220nm	125 mL bottle	3	Y	
20-300 nm Poly PSL (5E8/mL > 30 nm nominal)	KFMT cocktail 20-300nm	125 mL bottle	3	Y	
10 nm colloidal silica (1E11/mL)	KFMT CS 10nm	125 mL bottle	3	Y	diafiltered
20 nm colloidal silica (1E11/mL)	KFMT CS 20nm	125 mL bottle	6	N	diafiltered
30 nm colloidal silica (1E11/mL)	KFMT CS 30nm	125 mL bottle	6	N	diafiltered
50 nm colloidal silica (1E11/mL)	KFMT CS 50nm	125 mL bottle	6	N	diafiltered
100 nm colloidal silica (1E11/mL)	KFMT CS 100nm	125 mL bottle	6	N	diafiltered